

AB Max™



from  **Grafton Optical**

Advanced anterior blepharitis treatment

Feel the power of the pulse
with AB Max™



01923 233 980
sales@graftonoptical.com

Exclusively available in the UK through Grafton Optical


AB Max™ pulse mode provides the most effective and efficient treatment for anterior blepharitis on the market

AB Max™ is an advanced **5-8 minute deep cleansing treatment** to remove bacteria, crust and debris from the anterior eyelid margin and eyelashes. AB Max™ provides **advanced functionality** utilising two on-board computers and proprietary tips, specifically designed to treat Anterior Blepharitis. The AB Max™ **unique PULSE mode** provides improved patient outcomes and a more comfortable treatment experience.



To use, simply remove the AB Max™ handset from the charging base, open a new procedure pack and attach one of the two AB Max™ Tips. Fill the reservoir on the charging base with your liquid lid scrub solution and soak the AB Max™ Tip for 10-15 seconds.

Press the power button on the handset once for forward mode, again for reverse mode, and a third time for pulse mode. Clean the upper or lower outer lid margin for 20-30 seconds using each mode. Remove the AB Max™ Tip and attach a new one to treat the second eye.



AB Max™ provides better outcomes and higher profits

It is estimated that **Blepharitis accounts for 1 in every 20 eye problems** reported to GPs. This represents a huge opportunity for opticians looking to offer effective treatment plans to their patients, or who wish to enhance their current offering of ocular wipes, drops and masks.

The AB Max™ provides the same forward and reverse functionality as competitor devices, as well as a **unique PULSE mode**. Specifically engineered to remove even the worst crust and debris while massaging the anterior eyelid margins, this provides better patient outcomes and a more comfortable treatment experience.

The AB Max™ not only provides you with the most advanced treatment for Anterior Blepharitis, it also adds to your profit margin and provides a great return on investment. Treatments are usually repeated at regular intervals, depending on the severity of the case.


Maximise your treatments Maximise your profits

- ✓ Better treatment outcomes
- £ Affordable for patients
- ↑ Higher profits
- 🕒 Fast ROI



	AB Max™	Competitor device
Forward and reverse modes	●	●
Rechargeable handpiece with charging base	●	●
Patent pending pulse mode	●	
Affordable handpiece/ charging base	●	
Affordable procedure packs	●	
3 year standard warranty	●	

*Anterior blepharitis is blepharitis on the outer eyelid margin, including the eyelashes.



AB Max™ is a product of Myco Industries, Inc. and is exclusively available in the UK through Grafton Optical. AB Max™ is for use by medical professionals only.

Unit 7 River Park Industrial Estate, Billet Lane
Berkhamsted, Hertfordshire, HP4 1TB
sales@graftonoptical.com - 01923 233 980